

Impression taking at implant level

PF Ø 3.5 mm

PF Ø 4.0 mm



PF Ø 4.5 mm

PF Ø 5.0 mm



PF Ø 6.0 mm

Screws are not included with impression copings. The screw length may be determined according to the height of the impression coping and the tray technique used (open or closed, see page 13).



Impression coping for open tray technique, single use, retentive, screw-retained, aluminum, anodized

		Height					
		external Ø					
	cylindrical	10.0 mm	3.8 mm	4.3 mm	4.8 mm	5.3 mm	6.3 mm
			3.03.010	3.03.018	3.03.011	3.03.012	3.03.013
	conical	10.0 mm		5.0 mm	5.5 mm	6.0 mm	7.0 mm
				3.03.078	3.03.071	3.03.072	3.03.073



Impression coping for open tray technique, reusable, retentive, screw-retained, titanium

		Height					
		external Ø					
	cylindrical	16.0 mm	3.8 mm	4.3 mm	4.8 mm	5.3 mm	6.3 mm
			3.04.020	3.04.021	3.04.022	3.04.023	3.04.024
	conical	16.0 mm		5.0 mm	5.5 mm	6.0 mm	7.0 mm
				3.04.037	3.04.038	3.04.039	3.04.040


Impression coping for closed tray technique, single use, repositionable, screw-retained, titanium

		Height					
		external Ø					
	8.0 mm	short	3.9 mm	5.0 mm	5.5 mm	6.0 mm	7.0 mm
			3.03.074	3.03.079	3.03.075	3.03.076	3.03.077
	12.0 mm	long					
			3.03.014	3.03.019	3.03.015	3.03.016	3.03.017

Impression coping for closed tray technique, multiple use, repositionable, screw-retained, titanium

		Height					
		external Ø					
	8.0 mm	short	3.9 mm	5.0 mm	5.5 mm	6.0 mm	7.0 mm
			3.04.030	3.04.031	3.04.032	3.04.033	3.04.034
	12.0 mm	long					
			3.04.046	3.04.047	3.04.048	3.04.049	3.04.050

Positioning handle, stainless steel

		Length		Ø	
	17.0 mm	6.0 mm			
			3.03.522*		

* The positioning handle 3.03.522 may be used for all impression copings at implant level (except impression copings 16 mm of height).

PF Ø 3.5 mm

PF Ø 4.0 mm

PF Ø 4.5 mm


PF Ø 5.0 mm

PF Ø 6.0 mm

Screw for impression copings, reusable, 4-lobe, stainless steel

Length

12.5 mm	extra short	3.03.572					
16.5 mm	short	3.03.573					
21.5 mm	long	3.03.574					
27.5 mm	extra long	3.03.580					
13.2 mm	extra short					3.03.575	
17.2 mm	short					3.03.576	
22.2 mm	long					3.03.577	
28.2 mm	extra long					3.03.581	



Selection of appropriate screw length:

