

PRESS RELEASE

Grenchen, 12th March 2018

Greater added value with individualized abutments

With the new milled abutments for CAD/CAM type II, Thommen Medical has been able to boost the production capacity of its existing machine park and thus its in-house added value. This creates new opportunities for more entrepreneurial potential both for dental technicians and for dental laboratories. The milled abutments for CAD/CAM type II are available immediately in platforms 3.5-6.0mm.

The milled abutments for CAD/CAM type II are used with the PreFace® Abutment Holders to produce individualized titanium abutments. The milled abutments for cone and telescope construction and parallelized abutments for cemented bridge constructions can also be used. The new milled abutments are now available as a set that includes the abutment screw. The CAD libraries are available at www.thommenmedical.com. You can find detailed information on the new milled abutments for CAD/CAM type II in the relevant "Instructions For Use".

The milled abutments for CAD/CAM type II are machined in the Original Medentika® PreFace Abutment Holder in a suitable milling system. PreFace Abutment Holders can be ordered directly from the milling machine manufacturer. You can also obtain information on the use of the Medentika® PreFace Abutment Holder from your milling machine manufacturer.

Compatible milling machines:

- Datron D5®
- Dental Concept DC1/DC5®
- imes-icore®
- MB Maschinen Cobra Mill®
- Primacon PFM 24 mediMill®
- R+K®
- Röders RXD®
- VHF®
- Wissner Gamma 202®



Contact

Thommen Medical AG
Neckarsulmstrasse 28
CH-2540 Grenchen
Tel. +41 61 965 90 20
Fax +41 61 965 90 21
presse@thommenmedical.com
www.thommenmedical.com

The press release is also available online:

<http://www.thommenmedical.com/en/company/media/press-releases.html>

About Thommen Medical

Thommen Medical is a Swiss designer and manufacturer of a dental implant system known for incredible precision and impactful design. We are unique with an incredible heritage of over 30 years of Swiss precision, innovation, and functional design that has formed the basis for our products. Of high quality yet simple to use, our products deliver predictable, aesthetic results without compromise. Thommen Medical's headquarters and private manufacturing facility are in Grenchen, Switzerland, a



location that is the heart of Swiss precision manufacturing. Thommen Medical is an international company with subsidiaries and distributors throughout the world.