

PRESS RELEASE

Grenchen, 22.11.2022

Thommen Medical presents the next generation of SPI®VARIOmulti tertiary components

SPI®VARIOmulti abutments offer an optimal solution for screw-retained, temporary, or definitive long-span bridge restorations. The system has been optimized for conventional and digital workflow and features a full range of tertiary components. Now, the next, improved generation of SPI®VARIOmulti tertiary components is being introduced, which includes three major innovations.

Optimized for the fabrication of temporary reconstructions, the temporary copings now feature full-length macro-retentions for optimal support. The next generation of titanium impression copings has an integrated predetermined breaking point, which enables fast, manual shortening in addition to cutting. Another innovation concerns the protective cap, which now consists of only one part with an integrated thread, resulting in a simplified application.

The next generation of SPI®VARIOmulti tertiary parts – for fast and optimal care.

Contact

Thommen Medical AG

Neckarsulmstrasse 28

CH-2540 Grenchen

Tel. +41 61 965 90 20

Fax +41 61 965 90 21

info@thommenmedical.com

www.thommenmedical.com

Information on Thommen Medical:

Thommen Medical is proud that leading dentists choose to work with the Thommen Implant System, as only the best dental implant system is good enough for their patients. The dental implant system from Thommen Medical is easy to use and means that first class clinical results with an excellent esthetic result can be achieved. Also, after more than 35 years of clinical experience and working in close collaboration with international leading dentists, dental technicians and well-respected academic institutions in research and development, Thommen Medical is continuing towards its objective that every patient should have access to the best dental restoration solutions available. Thommen Medical manufactures the products in its own factory in Grenchen, Switzerland, and has subsidiaries in Switzerland, Germany and in the USA. Thommen Medical is represented by distribution partners in other important markets in Europe, Asia and the Middle East.

Thommen Medical AG | Headquarters

Neckarsulmstrasse 28, CH-2540 Grenchen, Tel. +41 61 965 90 20, Fax +41 61 965 90 21, info@thommenmedical.com
www.thommenmedical.com

Driven by science, not trends.